

Plate Analysis Script Instructions

Version 1.0

Enhanced Instructions & Cleaned Reasoning

1. Core Fixes & Enhancements

- Correct the Red/Blue swap: assign color channels properly.
- Add three numeric columns: Red, Green, Blue.
- After printing the well-level table, display:

Mean \pm 2·SD for R: ..., G: ..., B: ...

2. Command-Line Modes

--mode all Analyze all 96 wells; output full plate table and overall RGB statistics.

--mode left Analyze wells A1–H6 (48 wells); output left-side table and its RGB statistics.

--mode right Use left side as control, compute its Mean \pm 2·SD for each channel. On the right side (wells A7–H12), add a `WithinRange` column (TRUE/FALSE) based on whether each well's RGB triplet lies within control \pm 2·SD. Then output the right-side table, control stats, and "Number within range: {count}".

3. Additional Features

- **--avoidedges**: exclude outer-border wells from analysis.
- **--sd_range K**: use $K \times \sigma$ instead of fixed 2·SD.
- **--stats [mean,median,p25,p75,IQR]**: select summary statistics.

4. Image Handling

- Standard naming: `img_001.png`, `img_002.png`, ...
- `auto_crop.py`: crops a center circle (radius = 20 px) around each well.
- See inline comments and `docs/tutorial.mp4` for capture standards.

5. Output & Export

- **--output results.tsv**: save the well-level table as TSV.
- **--summary results_summary.txt**: save the statistics report.

6. Advanced & Future Extensions

- Batch-mode: process an entire directory of plate images.
- Classifier training: build color-vs-dilution models from known data.
- Dataset builder: generate CSV of dilution factors vs. RGB metrics.
- GitHub integration: set up CI tests, release tagging, and documentation site.

7. Usage Examples

- `python plate_analyze.py --mode all --stats mean,median --avoid.edges`
- `python plate_analyze.py --mode left --sd.range 1.5`
- `python plate_analyze.py --mode right --output test_results.tsv`

8. Notes

- Images must reside in `./images/` following the naming convention.
- Adjust crop radius in `auto_crop.py` for precise centering.
- Refer to `docs/tutorial.mp4` for a complete screen-recorded walkthrough.